

# PHASE ZERO SENSORS

## FAQs

1. How to purchase sensors?
  - CTI accepts all major credit cards or PO's.
2. I want to create a customized sensor; how do I get started?
  - Option 1: On CTI's website, click 'Discuss Your Project' and fill out the form
  - Option 2: call 800.706.0618
3. What type of connector do we suggest?
  - Electrodes are designed to fit in 2mm pitch connectors.
4. What if I need additional capabilities or services?
  - Click on the 'Discuss Your Project' on CTI's website, fill out the information and one of our team members will contact you soon.
5. Why should I choose a carbon sensor?
  - Carbon sensors provide a better bonding of deposited reagent due to its topography.
6. Why are the gold sensors plasma treated?
  - Due to gold's unique material characteristics it removes impurities, makes the surface more hydrophilic and improves electrochemical response.